

# C100 Cardiovascular Monitor

## Specification



### Standard Configuration:

12 Lead -ECG,RESP, Dual-Temp, Comen SPO2,NIBP

### Optional Configuration:

Nellcor/Masiom SPO2, 4-IBP,EtCO2 module, ICG,AG,C.O. , C100A

Telemetry monitor ,built-in WIFI module

### Safety Standards

ISO 13485:2003 approved, CE marking according to MDD93/42/EEC

### Physical Characteristics

Product Size:	305mm×168mm×308.4mm
Net Weight:	3.5kg (standard configuration)
Display:	
Type :	12" color TFT touch screen
Resolution:	800×600
waveforms:	up to 12 waveforms
External display	1 display through VGA

### Operation Environment

Temperature:	5-40℃
Humidity:	lower than 80%
Power requirement:	
AC Volatage:	AC100-250V, 50/60Hz
Current:	1.9A~0.8A
Battery:	
Battery Type:	Lithium-ion
Battery Capacity:	2600mAh 16.8V
Battery working time	2 hours
Battery Recharging Time:	Maximum 4 hours for charging;
Date storage	
Alarm Event Recall:	500 groups

NIBP Recall:	2000 groups
Trend Graph:	150 hours
Trend Table:	150 hours
Alarm Event Recall: (+ARR recall )	500 groups
ECG recall	7×24 hours of 12 lead ECG (8G SD card)
EDG diagnostic report	3000 groups
<b>Interface</b>	
Connectors:	USB interface, and PCL printer can be connected
	RJ45 Ethernet interface
	Nurse Call interface.
	AC interface
	Grounding pole
	DVI-D interface
	Multi-function interface to connect with sub-monitor and plug in slot
	IR module interface, it supports the following modules: CO2, dual IBP, AG, ICG,C.O.

### Recorder:

Trace:	Waveform: 6 channels for patient monitor , 12 channels for ECG
Type :	Built-in; thermal array
Speed:	Speed:25mm/s,50mm/s
Record width:	216 mm
Printer:	External printer supported

**ECG**

Lead Type	Cardio Tec™ 12-leads ECG Analysis, 5 Lead and 3 lead selectable
Lead selection:	3-lead: I; II; III 5-lead: I; II; III; aVR; aVL;aVF; V 12-Lead I; II; III; aVR; aVL;aVF; V1-V6.
Gain Selection:	2.5mm/mV ,5.0mm/mV,10mm/mV, 20mm/mV, auto
Sweep Speed:	6.25,12.5, 25, 50mm/s
Bandwidth :	MON Mode: 0.5-40Hz DIA Mode: 0.05-130Hz OPE Mode:1-20Hz
Defib. Protection:	ST Mode: 0.05-40HZ Withstand 4000VAC/50Hz voltage in isolation ,Again electrosurgical interference and defibrillation
Recovery time :	3s
ST SEGMENT detection:	-2.0mV - +2.0Mv
Arrhythmia Analysis :	Yes
pacemaker detection:	detectable
Alarm:	Yes, audible and visual alarm, alarm events recallable
12 lead ECG Analysis:	yes

**Heart Rate**

Range:	Adult: 15-300bpm Pediatric:15-350bpm Neonate:15-350bpm
Resolution:	1 bpm
Accuracy:	±1% or ±1bpm (whichever is greater)

**Respiration**

Method:	impedance method
RR measurement range:	Adult: 6-120bpm Pediatric/Neonate: 6 -150bpm
Resolution:	±1 bpm
Accuracy:	6-150rpm:±1rpm or 1%, whichever is greater 0-6rpm: unspecified
RESP Apnea Alarm:	Adult :10s-60s Child and Neonate: 10s-20s
Alarm:	Audible and visual alarm; alarm events recallable
Sweep Speed:	6.25,12.5,25 mm/s
Gain Selection:	X0.25, X0.5, X1, X2,X3,X4,X5

**Nellcor SPO2**

Measurement & alarm	0~100%
---------------------	--------

range:

Resolution:	1%
Accuracy:	±2% (70-100%, MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I and MAX-FAST sensors); ±3% (70-100%, D-YS, DS-100A, OXI-A/N and OXI-P/I sensors); 0-69% unspecified

PR Measurement Range: 20~300bpm

Resolution:	1bpm
Accuracy:	±3bpm
Alarm range:	20-250bpm

**Masimo SPO2**

Measurement & alarm range	1~100%
Resolution:	1%
Accuracy:	±2% (70~100%,Ped, non-motion+ ±3% (70-100%, Neo, non-motion); 0-69% unspecified
Alarm range	0~100%
Pulse rate:	Range: 25~240bpm Resolution:1bpm Accuracy: ±3bpm(non-motion) ±5bpm(motion) Alarm range:25~240bpm

**Comen SPO2**

Measurement & alarm range	1~100%
Resolution:	1%
Accuracy:	±2% (70~100%,Ped, non-motion+ ±3% (70-100%, Neo, non-motion); 0-69% unspecified
Alarm range	0~100%
Pulse rate:	Range: 20~300bpm Resolution:1bpm Accuracy: ±3bpm(non-motion) ±5bpm(motion) Alarm range:20~300bpm

**NIBP**

Method	Automatic oscillometric
Work mode:	Manual / Automatic/Continuous
Measurement Unit:	mmHg / Kpa selectable
Measurement types:	Systolic, Diastolic, Mean
Range of systolic pressure:	Adult Mode:40-270mmHg Pediatric Mode:40-200mmHg

Range of diastolic pressure:	Neonate Mode 40-135mmHg Adult Mode:10-215mmHg Pediatric Mode:10-150mmHg
Range of mean pressure:	Neonate Mode 10-100mmHg Adult Mode:20-235mmHg Pediatric Mode:20-165mmHg Neonate Mode 20-110mmHg
Over-pressure protection:	both Hardware and software over pressure protection
Accuracy:	Less than $\pm 5$ mmHg
Resolution	1mmHg
PR from NIBP:	40-240bpm

#### Temperature (Dual Channel)

Measurement range:	0-50°C (32 to 122F)
Resolution:	0.1°C
Accuracy:	$\pm 0.1$ °C (without probe)
Channel:	2 channels / T1, T2, TD (Temperature Difference)

#### IBP

Channel	4 channels
Measured Pressure	ART,PA,CVP,RAP,LAP,ICP,P1, P2
Measurement Range:	-10mmHg~300mmHg
Alarm Range:	-50mmHg~300mmHg
Accuracy:	$\pm 2\%$ or $\pm 1$ mmHg, whichever is greater

#### EtCO2(Sidestream)

CO2 Measurement Range:	0 -150 mm Hg, 0 to 79%, 0 to 20kPa (at 760mmHg)
Accuracy:	$\pm 2$ mm Hg (0 – 40 mm Hg) $\pm 5\%$ of reading (41 – 70 mm Hg) $\pm 8\%$ of reading (71 –100 mm Hg) $\pm 10\%$ of reading (101 –150 mm Hg)
awRR range	0-150rpm
RR accuracy:	$\pm 2$ rpm (0-70rpm) $\pm 5$ rpm (>70rpm)
Response time:	<240msec (10% to 90%)
Delay time:	<2s

#### EtCO2(mainstream)

CO2 Measurement Range:	0 -150 mm Hg, 0 to 79%, 0 to 20kPa (at 760mmHg)
Accuracy:	$\pm 2$ mm Hg (0 – 40 mm Hg) $\pm 5\%$ of reading (41 – 70 mm Hg) $\pm 8\%$ of reading (71 –100 mm Hg) $\pm 10\%$ of reading (101 –150 mm Hg)
awRR range	0-150rpm
awRR precision	$\pm 2$ rpm (0-70rpm) $\pm 5$ rpm (>70rpm)
Delay time:	<2s

#### Multi-Gas/O2

Method:	Infrared Absorption
Gas sorts:	CO2,O2,N2O,Des,Iso,Enf,Sev,Hal, (optional paramagnetic sensor)
Measurement range:	CO <sub>2</sub> : 0-15 vol% $\pm(0.2$ vol% + 2 % of the reading) N <sub>2</sub> O: 0-100 vol% $\pm(2$ vol% + 2 % of the reading) HAL, ISO, ENF: 0-8 vol% $\pm(0.15$ vol% + 5 % of the reading) SEV : 0-10 vol% $\pm(0.15$ vol% + 5 % of the reading) DES: 0-22 vol% $\pm(0.15$ vol% + 5 % of the reading) O <sub>2</sub> : 0-100vol% $\pm(1$ vol% + 2 % of the reading)
Data output:	Fi and Et values
Respiratory rate:	0-150rpm
Preheating time:	ISA CO <sub>2</sub> : < 10s, ISA OR+/AX+: < 20s
Alarm:	User selectable alarm limits for all measurement
others:	Up to 4 waveforms displayed Agent mixture detection MAC value displayed

#### Cardiac Output(C.O.)

Method:	Thermodilution
Measurement range:	C.O.: 0.1~20L/min BT: 25~43°C IT: 0~27°C
Resolution:	C.O.: 0.1L/min BT、IT:0.1°C
Accuracy:	C.O.: $\pm 5\%$ or $\pm 0.1$ L/min, whichever is greater BT、IT: $\pm 0.1$ °C (no sensor)
Alarm Range:	BT Hi limit: (LO limit +0.1) ~43°C (C.O.) BT Lo limit:23.0~(Hi limit-0.1) °C Step: 0.1°C

#### Impedance Cardiography (ICG)

Method:	Thoracic Bio-impedance measurement
Measurement range:	SV: 5~250 ml/b HR: 40~250 bpm C.O.: 1.4~15L/min
Accuracy:	SV: unspecified HR: $\pm 2$ bpm C.O: unspecified
Alarm range:	CI: HI (LO+0.1) ~15.0 L/min/m <sup>2</sup>

LO  $0 \sim (HI - 0.1) L/min/m^2$   
 TFC: HI (LO + 1)  $\sim 150/K\Omega$   
 LO  $10 \sim (HI - 1) K\Omega$

Alarm Deviation

CI:  $\pm 01. L/min/m^2$   
 TFC:  $\pm 1 K\Omega$

COMEN, one of the leading medical device manufacturers in China, and admitted to be the quickest developing company in medical industry.

COMEN design and produce medical devices for OR, NICU, and ICU. The product range covers Anesthesia devices, Monitoring & diagnostic products, Infant incubators, ventilators, medical pendants, operation light, all designed with specific patient data managing system to link-in to multiple medical database.

Based in Shenzhen - the hi-tech "headquarter" of China, COMEN has 30 branch offices operating business in China and 8 international branch offices in: USA, Russia, Turkey, Mexico, Thailand, India, Malaysia and Indonesia.

With all over 50,000 install bases in hospitals and medical centers globally, COMEN products are accepted and favored with its friendly design and excellent performance.

Shenzhen Comen Medical Instrument Co.,Ltd

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Gongming Sub-district, Gongming New District, Shenzhen, P.R.China

Tel: +86-755-2640 8879 2641 0713

Fax: +86-755-2643 1232

Website: en.comen.com

E-mail: info@szcomen.com

Specifications subject to change without notice

CFDA CE 0120

